

# A versatile, yeast-based synthetic heavy chain-only antibody platform that facilitates the construction and functional screening of multispecific CD28 T cell engagers

Ean L. Small, Alec Desai, Nathan Sharkey, Melissa Durkin, Cory Ahonen, Tushar Jain, Bianka Prinz, and Bradley Lunde



## BACKGROUND

Heavy chain-only antibodies (HCAbs) are emerging as key modalities in a variety of therapeutic contexts due to their unique features. Notably, since HCAbs do not require a light chain, they enable the efficient generation of bi- and multi-specific antibodies. We have generated a human, synthetic HCAb discovery platform in yeast, which has been utilized in efforts against multiple therapeutically relevant targets, including the isolation of CD28-binding HCAbs. The variable domains from our CD28 HCAbs have been subsequently engineered into IgG-like bispecific antibodies where one arm contains the CD28 single domain antibody, and the other arm contains a Fab against a tumor-associated antigen (TAA). These bispecific costimulatory molecules show enhancement of T-cell activation and cell killing when combined with CD3 x TAA bispecific molecules.

### Adimab synthetic HCAb library design and platform overview

Figure 1: Library design principles





Figure 2: HCAb discovery project flow



## RESULTS

### MACS/FACS-based selections using yeast HCAb platform delivers robust HCAb outputs for a variety of therapeutic targets

Figure 3: Binding affinity and epitope coverage

Platform delivers HCAbs with nanomolar binding affinity (A) and broad epitope coverage (B)






Figure 4: Affinity maturation yields large affinity improvements against Hen Egg Lysozyme (HEL) (A) and three undisclosed targets (B)




### Output HCAbs show favorable developability properties

Figure 5: Favorable developability for selected HCAbs across multiple Adimab campaigns




### Isolation of CD28-binding HCAbs

Figure 6: Overview of the CD28 HCAb selection process. HCAb libraries were subjected to MACS and FACS selections using CD28 antigen; isolated antibodies were characterized; a panel of leads were selected as candidates for multi-specific antibodies




### Characterization of CD28-binding HCAbs

Figure 7: Affinity & cell binding, developability, and functional activity of CD28 HCAb outputs



### Efficient, high-throughput generation of multi-specific panels


Figure 8: Adimab Chain Exchange (ChEx) is applied to generate panels of multi-specifics. ChEx generates panels of multi-specific antibodies for functional assays (A); analytical Ion Exchange (aIEX) output of CD28 x HER2 bispecific (B)



### CD28 x HER2 bispecifics show costimulatory activity

Figure 9: Adimab αCD28 in bispecifics

Experimental set-up (A); Adimab αCD28 bispecifics enhance SKOV-3 cell killing (B) and increase cytokine release (C) when paired with αCD3 panel



## Conclusions

- We have designed and deployed a novel, synthetic HCAb platform that allows for facile isolation of multiple nanomolar binders against diverse therapeutic targets
- Robust library design resulted in HCAb outputs with favorable developability properties: favorable polyspecificity and hydrophobicity profiles, high SEC purity, and a range of sdAb melting temperatures
- Robust design of affinity maturation libraries improves binding affinities while maintaining favorable developability properties
- Generated panel of CD28-binding HCAbs that show stimulation of T-cell activity via IL-2 secretion and promoter activation
- CD28 bispecific T-cell engagers show enhancement in T-cell activation when combined with our CD3 x TAA bispecific molecules

## Acknowledgements

We would like to thank the multiple teams within Adimab that contributed to this work including all members of High-throughput Expression, Protein Analytics, Computational Biology/Platform Technologies, Core, Mammalian, and Antibody Engineering for their contributions to and continued use of the platform.

